

# fx SEM™

**POWERFUL TECHNIQUE,  
QUICKLY ANALYZE,  
RESULTS IN SECONDS**





The new fX SEM™ custom x-ray source is designed exclusively for use on electron microscopes. The compact design, and slide mounting, allow very close coupling to the sample. The orientation yields high “flux” (x-rays) in small to large excitation areas on the sample surface. The fX SEM™ offers excitation areas 500μ to 25mm. The integrated high-voltage power supply operates up to a maximum power of 10 watts (35 kV and 0.1 mA depending on anode material). The close coupling provides XRF analytical results comparable to those from traditional “benchtop” or “standalone” units. The fX SEM™ is designed so that it does not interfere with the normal operation of the electron microscope, including the use of the electron beam on the same sample, at the same time collecting all elements simultaneously. No special cooling is required.

## Xb TECHNICAL SPECIFICATIONS


<b>Anode Type</b>	End-window transmission
<b>Target Material</b>	Ag, Mo & W
<b>Accelerating Voltage</b>	10-35kV
<b>Beam Current</b>	0-100μA
<b>Anode Spot Size</b>	<500μm
<b>Collimator Size</b>	200μm, 500μm, 1000μm (others available)
<b>Source Filters</b>	Available upon request
<b>Cooling Requirements</b>	Conduction cooled, no fan required
<b>Controls/Safety</b>	Variable control kV/μA, X-ray on/off buttons, kV/μA display. Interlocked to SEM, keyed power-on switch, Integrated high-voltage power supply, HV-On lamp, warning beacon

### But why?


Electron beams (from scanning electron microscopes) produce very high backgrounds hiding the trace elements in the sample. X-rays, from a true “x-ray” source don’t have this effect. Using the fX SEM™ low ppm levels of elements can be easily identified, quantified, and even producing trace level x-ray maps to view elemental distribution of trace elements in your sample.

#### INCREASE DETECTION

EDS Only




With SEM-XRF




#### IMPROVE ACCURACY

EDS Only




With SEM-XRF




#### SEPARATE PEAK OVERLAPS

EDS Only



With SEM-XRF



IXRF Systems, Inc.  
10421 Old Manchaca Rd., Suite 620  
Austin, TX 78748

Telephone: 512.386.6100  
Fax: 512.386.6105  
Email: info@ixrfsystems.com