



EDS SYSTEM UPGRADES

**BECAUSE YOUR EQUIPMENT
IS NEVER TOO OLD TO UPGRADE**



UPGRADES




Have a working detector and just need new software and electronics? No problem, IXRF replaces everything but the existing detector. Our electronics interface with every EDS detector. IXRF guarantees resolution and light element performance. Advanced Ethernet electronics enhance light element performance and allow faster digital imaging and x-ray elemental mapping. As a testament to our confidence in quality, IXRF includes a one (1) year warranty on our electronics.

SOFTWARE UPGRADES

IXRF customers never pay for software upgrades! It's our thanks for analyzing with us. New feature developments are made available, at no cost, for existing primary customers for the life of their system. We have carried our customers through Windows 3.11, 2000, NT, XP, Vista, and now Windows 7 at no additional charge to anyone in any country!

DETECTOR UPGRADES

IXRF can rebuild, maintain, and support, your existing detector. Si-Li crystals can be upgraded. Light element windows replace Be windows for light element detection. For customers that would like to move to SDD (Silicon Drift Diodes), we offer these as well.



FEATURES

- SPECTRA
- MAPPING
- LINESCANS
- FULL AUTOMATION
- ALL INCLUSIVE
- FREE UPGRADES
- AND MORE



IXRF Systems, Inc.
10421 Old Manchaca Rd., Suite 620
Austin, TX 78748

Telephone: 512.386.6100
Fax: 512.386.6105
Email: info@ixrfsystems.com